


Precision Balance CENT



Application

CENT® Balances provides maximum weighing performance and full network capability, they have excellent weighing performance.

Features

- Large easy to read LCD display
- RS232 interface for connecting computers, printers and other equipment to facilitate the recording and processing of data
- High Precision Loadcell
- Optional battery power
- Calibration weight is included

TECHNICAL SPECIFICATION

Model	CENT600	CENT1000	CENT2000	CENT3000
Weighing Capacity:	600g	1000g	2000g	3000g
Readability:	0,01g			
Repeatability:	± 0,2mg			
Linearity	± 0,2mg			
Stabilisation time:	3-5 seconds			
Calibration:	External			
RS232 Interface:	Yes			
Mode:	Weighing; parts counting; percentage;			
Dimensions (L * W * H):	340 * 215 * 350			