

Analytical Balance ANBEC



Application

Analytical Balance External Calibration is ideal for laboratory and general purpose weighing.

Features

- Large easy to read LCD display
- overload protection
- automatic fault detection
- unit conversion
- RS232 interface for connecting computers, printers and other equipment to facilitate the recording and processing of data
- Easy to use, sealed keypad

TECHNICAL SPECIFICATION

Model	ANBEC100	ANBEC160	ANBEC210	ANBEC220
Weighing Capacity:	100g	160g	210g	220g
Readability:	0,1mg			
Repeatability:	± 0,2mg			
Stabilisation time:	3-5 seconds			
Calibration:	External			
RS232 Interface:	Yes			
Mode:	Weighing; parts counting; percentage;			
Dimensions (L * W * H):	340mm x 215mm x 350mm			