



HYBADEC Hydrostatic - Density Balance SET

External Calibration

€

- Capacity 2/3/5 kg
- Resolution 0,1g/0,01g
- LCD Display with backlight

Download
Catalogue



HYBADEC

Hydrostatic Density Balance



External Calibration

External Linear Calibration



High Precision Load Cell

High measuring range up to 5kg is provided with high capacity load cell,



Software for data collection

Optional software for communication with the scale provides the ability to collect data, control the measurement process ...

5000g/0,01g

**External
Calibration**

Big Pan Size

**Hydrostatic
SET**

**High Precision
Load Cell**

**Hydrostatic
Balance**