

ROVA901 Laboratory Rotary evaporator

Product Details

Product Number

PRO3502



Technical Specification:

Technical Specifications

Technical Specifications

Cooling type	Vertical
Cooling surface	1500 cm ²
Motor principle	DC motor
Speed range	5–300 rpm
Speed tolerance	Set rotation speed > 100 rpm ± 1%
Speed tolerance	Set rotation speed < 100 rpm ± 1 rpm
Lift type	Manual
Stroke	120 mm
Heating temperature range	Ambient to 180 °C
Temperature set resolution	1 K
Max filling volume	4 L
Heating power	1400 W
Input power	1500 W

Dimensions (W×D×H)	510 × 490 × 345 mm
Weight	25.76 kg
Ambient operating temperature	5–40 °C
Permissible relative humidity	80%
Protection class	DIN EN 60529 IP20
Power supply voltage	100–240 V
Frequency	50 / 60 Hz
DC voltage	24 V DC
Documentation	Catalog (PDF)

Accessories: