



LACE-12CH Hematocrit Centrifuge



LED Display




Neutral Temperature

Download Catalogue

Hematocrit Rotor

LACE12CH Hematocrit centrifuge with 24 capillary rotor included



Max. Capacity 4x180ml

Swing out Rotor



12000rpm/15800xg

Maximum speed and rcf meet all demand for clinical centrifuge



24 Capillary Rotor included



LED Display

Speed/RCF, Time, Parameters can be set during the centrifuge work..



Brushless Motor

COLO brushless motor is maintenance free, no carbon dust polution



Safety Function

Electric Lead Lock with alarm with safety function, disbalance detection, over heating of centrifuge protection

PRODUCT CODE	LACE-12CH
Max. Speed	12000rpm/15800xg
Max. Capacity	24x Hematocrit
Timer	0-99min
Noise	<45dBA
Display	LCD – Display: Speed, RCF, Time
Power supply	220V (±10%) 50-60Hz
Safety	Electro lid lock
Weight NW/GW kg	12/15kg
Product Dimension	355×270×205mm