



LACE-5RH Refrigerated Low Speed High Capacity Centrifuge

Angle/Swing out Rotor



LED Display



Refrigerated
-20°C to 40°C

Download
Catalogue



LACE32 is an economical universal centrifuge suitable for large clinical laboratories.

31392 Swing out Rotor 4 x 300ml 4000rpm x 2920xg

31382 Swing out Rotor 32 x 15ml 4000rpm x 3200xg

31383 Swing out Rotor 48 x 10ml 4000rpm x 2880xg

30690 Swing out Rotor 72 x 10ml 4000rpm x 3140xg

30589 Swing out Rotor 88x5ml 4000rpm x 3140/2840xg

31491 Micro Plate Rotor 2 x 4 x 96 wells 4000rpm x 2490xg

30641 Angle Rotor 12 x 50ml / 24 x 50ml 5000rpm x 3860xg

30614 Angle Rotor 6 x 100ml / 4 x 100ml 5000rpm x 3130xg



Max. Capacity 4x300ml

Swing out Rotor



5000rpm/4730xg

Maximum speed and rcf meet all demand for clinical centrifuge



Swing out/Angle rotors

Adapters available for this rotor



LED Display

Speed/RCF, Time, Parameters can be set during the centrifuge work..



Brushless Motor

COLO brushless motor is maintenance free, no carbon dust pollution



Safety Function

Electric Lead Lock with alarm with safety function, disbalance detection, over heating of centrifuge protection

Item Name	LACE-5RH
Max. Speed (Swing out rotor)	5000rpm/4730xg
Max. Capacity (Swing out rotor)	4 x 300ml
Noise	<58dBA
Timer	99H, 59min
Rotor type	Swing out and Angle Rotor
Power supply	AC 220V/50Hz
Product dimension	650 × 665 × 400mm
Weight kg NV/GW	118kg/130kg