

TITBUR2000 Electronic Burette



Download
Catalogue



High quality motor with excellent precision and accuracy



Electronic control provides fatigue-free operation



Remote control panel prevents manual disturbance during operation



Built-in magnetic stirrer for a complete and easy operation



Large LCD display parameters

- Digital Burette – TITBUR2000 provides accurate, precise and convenient bottle-top titration with optimum operator safety. It combines both magnetic stirrer and Burette functions , making it easier for the user to achieve reproducible results in the field of analytical chemistry, food industry, water analysis...

Features

- High quality motor with excellent precision and accuracy
- Electronic control provides fatigue-free operation
- Remote control panel prevents manual disturbance during operation
- Built-in magnetic stirrer for a complete and easy operation
- Large LCD display parameters

Type	TITBUR2000/50	Titbur2000/100
Max volume ml	50	100
Resolution	10µL	
Volume Accuracy	R=0.2% CV=0.07%	
Measurement error	according to ISO 8655-3	
Speed	16 Levels	
Operating Temperature Range	10°C-30°C	
Control Type	External control, including the stirrer and the burette program	

COLO Laboratory Equipment
 Polje ob Sotli 4
 SI-3255 Slovenia
 Tel.+386 64 222 724
 Email: sales@colo.si
 www.colo.si