



Ice maker

LFIM Laboratory Flake Ice Maker

Download
Catalogue



LFIM



High-quality stainless steel (housing and storage tank)
corrosion and durable



Unique water tank float system to ensure no residual water
over water without de-icing process



The small ice particles form the amorphous ice snow, ice forms small, can
penetrate narrow gaps fast cooling effect of ice bath, designed for laboratory
design

Type	Ice volume (Kg/24h)	Storage capacity	Ice shape	Dimension cm L*W*H	Input power W	NW-Kg	Refrigerant
LFIM20	20	10	Flakes	295×490×670	240	32	R134a/80g
LFIM30	30	10		350*525*670	240	35	R134a/80g
LFIM40	40	10		380×543×722	240	40	R134a/80g
LFIM50	50	10		380×544×723	240	45	R134a/80g
LFIM70	70	20		380×544×723	420	50	R404a/120g
LFIM85	85	20		548×611×1050	420	57	R404a/120g
LFIM100	100	25		548×611×1050	600	60	R404a/160g
LFIM150	150	35		548×611×1050	700	65	R404a/160g
LFIM200	200	55		680×641×1100	800	88	R404a/220g
LFIM300	300	55		720×641×1150	1300	100	R404a/220g
Power supply	220V +/-10% 50/60Hz						

Please note
 If this Ice maker does not match the specific needs of your application, or some options are not listed for sale, please feel free to contact us. Our manufacturing engineers will come up with technical solutions to meet your needs. We reserve the right to change technical specifications at any time.

COLO LabExperts 

Laboratory Equipment production and Distribution
 www.colo.si
 Polje ob Sotli 4
 SI-3255 Slovenia
 tel. +386 590 48 880 E-mail: office@colo.si