



UCB5 Ultrasonic Cleaning Bath

Price - €

Download
Catalogue



Digital Display

Two independent heating plate with analog regulation ...



Timer

Timing range 1-600 min adjustable



Stainless Steel Basket

Designed especially for ultrasonic cleaner to reduce ultrasound cover



Stainless Steel Structure

The shell and cover, made from stainless steel isolated to decrease ultrasound noises



Heating

Model	Inner Tank	Volume	Frequency	Power	Adjustable time	Drainage	Net weight	Packing size	Rack
	L/W/H (mm)	L	KHz	W	min		Kg	mm	
UCB-50	150*85*65	1	40	50	1-120	/		220*150*210	√
UCB-80	155*140*100	2	40	80	1-120	/	2.2	240*230*300	√
UCB-100D	300*150*100	4.5	40	100	1-600	/	4.5	400*240*435	√
UCB-120D	240*140*100	3	40	120	1-600	/	4.5	325*230*300	√
UCB-120E	240*140*150	5	40	120	1-600	/	4	325*230*350	√
UCB-3200D	300*155*150	6	40	180	1-600	/	5	400*240*380	√
UCB-5200D	300*240v150	10	40	200	1-600	√	6.8	450*325*375	√
UCB-5200E	300*240*150	10	40	250	1-600	√	7	450*325*375	√
UCB-5200F	300*240*150	10	40	300	1-600	√	7.2	450*325*375	√
UCB-4200D	300*240*200	14.4	40	400	1-600	√	9	445*320*460	√
UCB25-12D	500*300*150	22.5	40	500	1-600	√	13	640*395*375	√
UCB25-12E	500*300*150	22.5	40	600	1-600	√	13.5	640*395*375	√