

# MANIVAC3S Vacuum Membrane Filtration Manifold 3 place 300ml

## Product Details

Product Number

PRO5023



## Technical Specification:

- Prepare and detect and quantify microorganisms in liquid samples
- Membrane filtration
- 3 branch funnels
- Spin lock connection provide fast assembling and dissembling
- Accept 47mm and 50mm membrane filter
- Filter funnel is detachable, easy to sterilize by flame, autoclave or drying oven
- Funnel capacity 300ml
- Body, control valve, filter funnel, funnel support base: stainless steel SUS316

## Accessories: